

CAPSULES

preliminary notes and applications from Bioanalytical Systems, Inc.

Determination Of Oxytocin In Pharmaceutical Samples

Purpose

Quantitative analysis of oxytocin (F1) in pharmaceutical samples. A rapid and accurate method of determining the amount and purity of oxytocin synthesized by the pharmaceutical industry was needed.

Existing Methods

Bioassays based on administering oxytocin compounds to animals and measuring some aspect of their physiological effect, such as blood pressure, uterine contractions, or milk ejection pressure. These are tedious, imprecise, and time-consuming.

Reference

Quantitative Analysis of Oxytocin in Pharmaceutical Preparations by High-Performance Liquid Chromatography, M. Ohta, H. Fukuda, T. Kimura and A. Tanaka, J. Chromatogr., 402(1987) 392-395.

Conditions

Detector: UV, 210 nm

Sensitivity: 0.04 a.u.f.s.

Column: 5 μ m Zorbax TMS (250 x 4.6 mm),
temperature controlled at 40°C

Mobile Phase: acetonitrile: 50 mM phosphate buffer
(pH 5.0) 18:82. Flow rate was 1 mL/min.

Linear Range: 0.2 - 20 I.U./mL

Sample Preparation

Commercially prepared oxytocin and ethyl p-hydroxybenzoate (internal standard) were used. Injection volume was 80 μ L.

Notes

Measurement of oxytocin in commercial preparations by LCUV (F2) and bioassay (chicken blood pressure) produced nearly identical results.

The determination of oxytocin presented in this report can be duplicated using the BAS 400 Liquid Chromatograph.

The information in this publication is supplied as a service to our customers. Performance of the methodology has not necessarily been verified by BAS technical staff.




Figure 1. Structure of Oxytocin.




Figure 2. Chromatogram of (1) oxytocin, (2) chlorobutanol preservative, and (3) ethyl p-hydrobenzoate.

COPYRIGHT 1987, Bioanalytical Sysstem, Inc.