

# CAPSULES

preliminary notes and applications from Bioanalytical Systems, Inc.

## Phenobarbital in Serum

### Purpose

Determination of phenobarbital in serum.




Figure 1. Structure of phenobarbital.

Phenobarbital (F1, 5-ethyl-5-phenyl-2,4,6(1H,3H,5H)-pyrimidinetrione) is a widely used barbiturate classified as an anticonvulsant and sedative-hypnotic. It is used in the treatment of seizures (including those of the newborn) and as a long-acting sedative. Therapeutic doses range from 10-40  $\mu\text{g/mL}$  blood. Side effects associated with the use of phenobarbital may include impaired coordination (ataxia), liver damage, and drowsiness.

### Existing Methods

GC coupled with chemical ionization mass spectrometry, enzyme immunoassay and LC. Numerous methods have been published and are available commercially.

### Conditions

Detector: BAS UV-108 variable wavelength UV detector (220 nm)

Column: 3  $\mu\text{m}$ , C 18 reverse-phase, 100 x 3.2 mm (PN MF-6213)

Mobile Phase: 80% (v:v) 25 mM  $\text{Na}_2\text{HPO}_4$ , pH 7.0, 20% acetonitrile. Flow rate was 1 mL/min.

Detection Limit: 175 pg injected standard, 35 ng/mL serum (S/N = 3).

Linear Range: 0.5-100 ng injected standards, 5-40  $\mu\text{g/mL}$  serum.




Figure 2. Calibration curve for spiked serum samples.


### Sample Preparation

1. Prepare Bond-Elut® C 18 solid-phase extraction columns by washing with 1 mL methanol followed by 1 mL water.

2. Load the following onto each column: 50  $\mu\text{L}$  serum, standards as appropriate, and 0.05 M  $\text{NaH}_2\text{PO}_4$  (pH 5.0) for a final volume of 1 mL.

3. Wash the column with an additional 1 mL of phosphate buffer followed by two 1-mL washes of water.

4. Elute the samples with 1 mL methanol. Dry each in a stream of nitrogen or in a vacuum evaporator. Redissolve in 0.350 mL mobile phase and inject in 20  $\mu\text{L}$  aliquots.



**Figure 3.** Sample chromatograms of blank (left) and spiked (5 µg phenobarbital/mL) serum samples.

#### Notes

A calibration curve for spiked serum samples is presented in F2, and sample chromatograms are shown in F3.

Recovery of phenobarbital from spiked serum samples was 98%.

5-methyl-5-phenylhydantoin [1], guaifenesin [2] and 5-(p- tolyl)-5-phenylhydantoin [3] have been used as internal standards.

#### References

1. Bhargava, V.O., *J. Chromatogr.* 419 (1987) 421-425.
2. Shabbir, U. and F. Castro, *J. Liq. Chromatogr.* 9 (1986) 2269-2279.
3. Soto-Otero, R., E. Mendez-Alvarez and G. Sierra-Marcuno, *J. Liq. Chromatogr.* 8 (1985) 753-763.

