

# CAPSULES

preliminary notes and applications from Bioanalytical Systems, Inc.

## In Vivo Microdialysis Rat Brain - Purines

### Purpose

To determine free purines in the brain of freely-moving rats.

### Conditions

System: BAS-400 Purine Analyzer with  
Microdialysis System II

Detector: BAS UV-8, 254 nm

Column: Biophase Octyl, 5  $\mu$ m, 250 x 4.6 mm  
(BAS, PN MF-6032)

Temperature: Ambient

Mobile Phase: 10 mM NaH<sub>2</sub>PO<sub>4</sub> (pH 4.5) plus 9%  
CH<sub>3</sub>OH

Flow Rate: 1.2 mL/min

Detection Limit: 0.1, 0.2, 0.3, 0.4, and 0.6 ng,  
respectively, of hypoxanthine, inosine,  
guanosine, adenine, and adenosine at a S/N of  
3.

Microdialysis Probe: CMA/10, 3 mm membrane

Perfusion Speed: 2  $\mu$ L/min using CMA/100 Microin-  
jection Pump

Containment System: CMA/120 for awake rats


### Sample Preparation

Stereotactically implant a guide cannula into the rat caudate and cement it to the skull(1). Allow the rat to recover from surgery for several (2-3) days. The microdialysis probe is inserted while the rat is awake. Continuously perfuse Ringer's solution into the probe and collect dialysates in the CMA/140 fraction collector. Directly inject 5  $\mu$ L of each fraction into the Purine Analyzer.

### Notes

Chromatograms of purine standards and rat brain dialysates are presented in F1.

The determination of purines presented in this report can be duplicated on the BAS 200 Problem Solver.



**Figure 1.** Chromatogram of (A) 1  $\mu$ M purine standards, 10  $\mu$ L was injected, AUfs 0.02 (B) awake rat brain dialysates, 5  $\mu$ L was injected, AUfs 0.005 Peaks: 1. hypoxanthine 2. inosine 3. guanosine 4. adenine 5. adenosine

The information in this publication is supplied as a service to our customers.

### Reference

1. *Carnegie Medicin, Microdialysis User's Guide, 4th Edition, Stockholm, Sweden, 1988.*

COPYRIGHT 1988, Bioanalytical Systems, Inc.