

Detection of Amino Acids in Elastin

F3 shows the detection of amino acids in elastin. The gradient was modified to detect lysine and desmosine. These compounds are multiderivatized, and elute about 15 minutes after the last peak if the original chromatographic conditions are employed. The amino acid composition of elastin as determined by this method agreed well with those reported in the literature.




Figure 3. Detection of desmosine and isodesmosine in the elastin hydrolysate using gradient elution LCEC. Mobile phase A was 85% 0.005 M sodium citrate (pH 7.5), 10% methanol, and 5% THF. Mobile phase B was 10% 0.005 M sodium citrate (pH 7.5) and 90% methanol. Both mobile phase A and B were made 0.05 M with respect to sodium perchlorate. Gradient conditions were 50 to 100% B over 30 min; other conditions as in F2. Concentration of elastin, 4.0 $\mu\text{g/mL}$.

In addition, it has been found that the oxidation potential of the CBI-derivatized amino acids is strongly dependent on the R-group. Multiderivatized compounds, such as lysine, were found to have a lower oxidation potential than monoderivatized amino acids such as alanine. Therefore, voltammetry can be used for the identification of sample components. In reference 2, a dual electrode detector was used to compare the redox potential of the peaks coeluting with lysine and desmosine with that of derivatized standards. The voltammetry was found to be identical.

Notes

The derivatization volumes can be scaled down by a factor of 100. A large volume was used in this case because the sample quantity was not limited.

When determining trace levels of amino acids, it is important to thoroughly clean glassware and to use water that is amine free.

We have found that the addition of NaOH to the sodium borate increases the number of extraneous peaks. Therefore a 0.01-0.05 M sodium borate solution (unadjusted pH about 9.3) is used for all amino acid derivatizations.

It is very important to add the NDA last and to keep the sample protected from light.

References

1. P. de Montigny, J.F. Stobaugh, R.S. Glvens, R.G. Carlson, K. Srinivasachar, L.A. Sternson and T. Higuchi, *Anal. Chem.* 59 (1987) 1096-1101.
2. S.M. Lunte, T. Mohabbat, O.S. Wong, and T. Kuwana, *Anal. Biochem.* 178 (1989) 202-207.
3. M. D. Oates and J. W. Jorgenson, *Anal. Chem.* 62 (1990) 577-1580.

* U.S. PATENT # 4837166

BAS

