



CAPSULES

notes and applications from Bioanalytical Systems, Inc.

Determination of Purines in Microdialysates Using SepStik Columns

Purpose:

Determination of adenosine, guanosine and adenine (F1) in microdialysates.




Figure 1. Structures of adenosine (A), guanosine (B) and adenine (C).

In order to separate and sensitively detect purines in microdialysates, a SepStik microbore column was used. The 1 mm internal diameter increases the concentrations of the eluting purines up to 21-fold compared to standard LC columns.

Existing Methods

LCEC and LCUV with conventional columns.

Conditions

System: BAS 481 Liquid Chromatograph

Column: SepStik kit (BAS P/N MF-8949).

The packing was ODS 3 μ M silica in a 100 x 1.0 mm bed.

Mobile Phase: 7 mM NaH_2PO_4 containing 3.5% CH_3OH . Adjust pH to 3.04 after adding CH_3OH .

Flow Rate: 80 μ L/min actual measured flow rate.

For conventional pumps, use BAS P/N MF-8947 Flow Splitter Accessory kit, and maintain about 2 mL/min flow rate at the pump. Adjust splitter ratio to give 75-85 μ L/min through the microbore column by adjusting restrictor and/or pump flow rate.

Detector: BAS UV-108, with microcell.

Wavelength: 260 nm

Range: 0.002 AUFS

Rise Time: 1.0 sec.

Detection Limit: 50 pg injected yielded a S/N of 3.

The injection volume was 5 μ L.

Sample preparation:

Dialysate was directly injected onto the system.




Figure 2. Detection limit test

A: 500pg of each purine standard

B: Blank

C: 50pg of each purine standard

Peak Identification:

1. hypoxanthine
2. adenine
3. inosine
4. guanosine
5. adenosine




Figure 3. Purines in various rat striatum dialysates.

A: Blank (Ringer's solution)

B, C & D: Various dialysate samples

Peak Identification as in F2.