



## Determination of Salsolinol by Microbore LCEC

### Purpose

Develop a sensitive LCEC method for the determination of salsolinol (SAL) (**F1**) at a low concentration.




Figure 1. Structure of SAL

SAL was found to occur in human brain (1), and to be a precursor of 6, 7-dihydroxy-1, 2-dimethyl-1, 2, 3, 4-tetrahydroisoquinol (N-methyl-salsolinol), an endogenous neurotoxin (2). Therefore, a simple and reliable analytical method for the determination of SAL is necessary for studying the pharmacological effects of SAL.

### Existing Methods

Various analytical methods are available for the determination of SAL, such as, LCEC (3, 4) and GC-MS (5, 6). However, we needed to develop an LCEC method with a low detection limit. The following microbore LCEC method is relatively simple, rapid, and highly sensitive.

### LCEC Conditions

**System:** BAS 480e system equipped with an epsilon EC detector

**Electrode:** glassy carbon (3 mm)

**Potential:** +700 mV vs. Ag/AgCl

**Column:** UniJet, ODS, 3 $\mu$ m, 100 x 1 mm (BAS, P/N MF-8949)

**Mobile Phase:** 1L of buffer (0.1M monochloroacetic acid, 0.05 mM EDTA, 0.65 mM octylsulfate, sodium salt), pH to 3.1; 50 mL Methanol.

**Flow Rate:** 100  $\mu$ L/min

**Amount Injected:** 10  $\mu$ L

**Temperature:** 35°C

### Note

A representative chromatogram of SAL standard is shown in **F2**.




Figure 2. Typical chromatogram of 10 pg SAL

The calibration curves of SAL at the ranges of 0-200 pg and 0-1000 pg are presented in **F3a** and **F3b**, respectively. The detection limit is 2 pg on column.




Figure 3a. Calibration curve of SAL (0 - 200pg)

(continued)



## References

1. Sjoquist B., Eriksson A., and Winblad B., *Prog. Clin. Biol. Res.*, 90 (1982) 57.
2. Naoi M., Maruyama W., Dostert P., Hashizume Y., Nakahara D., Takahashi T., and Ota M., *Brain Res.*, 709 (1996) 285.
3. Baum S.S., Hill R., Kiianmaa K., and Rommelspacher H., *Alcohol*, 18 (1999) 165-169.
4. Deng Y., Maruyama W., Kawai M., Dostert P., Yamamura H., Takahashi T., and Naoi M., *J. Chromatogr.*, 689 (1997) 313-320.
5. Haber H., Roske I., Rottmann M., Georgi M., and Melzig M.F., *Life Sci.*, 60 (1997) 79-89.




Figure 3b. Calibration Curve of SAL (0 - 1000pg)