




Determination of Naringenin in Rat Plasma with the Culex® Automated Blood Sampler Coupled with LCEC

Purpose

Develop an automated blood sampling LCEC method for pharmacokinetic studies of Naringenin (NG-E, **F1**) in rat plasma.

NG-E is one of the active components in grapefruit juice that changes the pharmacokinetics of a variety of drugs. It also has the ability to activate polymorphonuclear (PMN) leucocytes to exhibit cytotoxic activity against tumor cells *in vitro*.



F1. Structures of NG-E.

Existing Methods

LCUV, LCEC, LCMS, GC/MS and capillary electrophoresis have been used to determine NG-E. Several of these methods suffer from inadequate lower limits of quantitation (LLOQ) and require tedious sample preparation steps. Our goal was a method suited to small blood volumes from rodent pharmacokinetic experiments.

Conditions

LCEC System: BASi 200B temperature control and solvent delivery system equipped with a LC-4C detector (BASi) for microbore chromatography

Electrode: glassy carbon working electrode, 3 mm (MF-1000)

Potential: +900 mV vs. Ag/AgCl

Guard Column: C18, 5 μ M, 14 x 1.0 mm (MF-8945)

Analytical Column: C18, 5 μ M, 150 x 1.0 mm (MF-8912)

Mobile Phase: 1 L of buffer containing 20 mM NaOAc and 27 μ M EDTA (pH 6.2) plus 300 mL of CH₃CN and 50 mL of CH₃OH

Flow Rate: 100 μ L/min

Injected Volume: 10 μ L

Temperature: 50°C

Detection Limit: 0.5 pg of NG-E on column

Quantitation Limit: 1 ng/mL of NG-E in rat plasma at S/N of 5:1

Blood Collecting System: freely moving rat containment device (Raturn®, BASi), automated blood sampler (Culex®, BASi) and refrigerated fraction collector (HoneyComb™, BASi).

Sample Preparation

Heparinized rat blood samples were centrifuged at 2000 g for 10 min. Each 60 μ L of plasma was transferred to a 1.7 mL centrifuge tube and extracted with ethyl acetate three times using 400 μ L, 240 μ L and 240 μ L, respectively. The combined organic solvent was dried with nitrogen gas at room temperature and the residue dissolved in 60 μ L mobile phase. Samples were diluted as needed to fall within the determined linear range. A 10 μ L aliquot was injected into the LCEC system.

Primary Animal Study

For the automated blood sampling experiment, a male Sprague-Dawley (250g) was implanted with a jugular vein catheter (CX-2010, BASi) and a femoral vein catheter (CX-2020, BASi). After surgery, the rat was installed in the Raturn, and allowed to recover for two days with free access to food and water. The rat was dosed at 5mg/kg of NG-E through the femoral vein catheter at 2mg/mL of NG-E in 30% cyclodextrin. Successive 120 μ L blood samples were collected through the jugular vein catheter into heparinized vials and kept refrigerated in a fraction collector according to a preset schedule.

Note

Chromatograms of NG-E standard (A), blank rat plasma (B) and plasma spiked with 10 ng/mL of NG-E (C) are shown in **F2**.

(continued on back)




The recovery of NG-E spiked into rat plasma at a concentration of 10 ng/mL, 50 ng/mL, 100 ng/mL, and 1000 ng/mL (n=5) was determined to be 77.11%, 74.50, 74.43%, and 82.40%, respectively.

Rat plasma spiked with NG-E was linear in the range of 25-200 ng/mL. The regression equation is $y=248386x-213486$ with $R^2=0.998$.


The developed method was used for a pharmacokinetic study of NG-E in a rat. **F3** illustrates the results of a single intravenous dose of NG-E administered to a rat.

References

1. H.W. Peng, F.C. Cheng, Y.T. Huang, C.F. Chen and T.H. Tsai., *J. Chromatogr. B*, 14 (1998) 369-374.
2. K. Ishii, T. Furuta and Y. Kasuya, *J. Agric. Food Chem.*, 48 (2000) 56-59.
3. S.H. Lee, Y.B. Park, K.H. Bae, S.H. Bok, Y.K. Kwon, E.S. Lee and M.S. Choi., *Ann. Nutr. Metab.*, 43 (1999) 173-180.
4. J.D. Spence. *Clin. Pharmacol. Ther.*, 61 (1997) 395-400.
5. T.B. Ng, F. Liu and Z.T. Wang., *Life Sci.*, 2000, 66(8):709-723.



F2. Chromatograms of standard NG-E (10 ng/mL) (A), extracted blank rat plasma (B), and extracted plasma spiked with 10 ng/mL of NG-E (C).



F3. Rat plasma NG-E concentration vs. time profile in a rat following a single intravenous administration of NG-E.